



Product Specification

Blood Gas-Monovette® , calcium-balanced lithium heparin, 1 ml, cap white/orange, connection: Luer (m)

Page 1



Product description

Order number	05.1146
Product description	Blood Gas-Monovette® , preparation: calcium-balanced lithium heparin, nominal volume: 1 ml, cap white/orange, label/print: transparent/orange, connection: Luer (m), (LxØ) without cap: 66 x 11 mm, 50 piece(s)/case, sterile, 50 piece(s)/case

Product characteristics

Type of collection	Arterial and venous
Type of preparation	calcium-balanced lithium heparin
Preparation concentration	50 IU/ml blood
Label/ Print	with plastic label
Colour of print/label	transparent/orange
Graduation	fill mark
Application	blood gas analysis
Connection	Luer (m)

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2022-11-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

Blood Gas-Monovette® , calcium-balanced lithium heparin, 1 ml, cap white/orange, connection: Luer (m)

Page 2

Size

Sample volume	1 ml
Diameter	11 mm
Length including cap	82 mm
Length excluding cap	66 mm
Length including cap and plunger	100 mm

Material & colours

Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	white/orange
Tube material	Polypropylene (PP)
Piston material	High Density Polyethylene (HD-PE)
Piston plunger material	Polystyrene (PS)

Purity & certification

Product category	In-vitro diagnostic, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

Packaging

Minimum order qty.	500
Type of inner packaging	case
Piece(s) / inner box	50
Piece(s) / case	500
Piece(s) / pallet	56000
Depth of box	145 mm
Width of box	108 mm
Height of box	63 mm
Depth of case	333 mm
Width of case	303 mm
Height of case	125 mm
Case volume	0.0126 cbm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2022-11-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

Blood Gas-Monovette® , calcium-balanced lithium heparin, 1 ml, cap white/orange, connection: Luer (m)

Page 3

Weight of product	0.0043 kg
Weight of case	2.69 kg
EAN of inner box	4038917247221
EAN case	4038917247214

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2022-11-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT